

PR515
Introducing the next generation of wellness positioning! The new PR515 is available in two sizes and comes standard with our PowerPillow™ and power lumbar support for even more comfort and positioning options. This luxurious power lift recliner features two incredible patents combining our extremely popular MaxiComfort Cloud™ positioning and Twilight Technology that provide superior comfort and quality found only at Golden. Achieve extreme Zero Gravity and TV positions while enjoying the world's finest rejuvenating seating experience!





PR-515-SME FABRICS

Calypso

Copper

Hazelnut

Shiraz

Sterling

Calypso

Hazelnut

Shiraz

Brisa Coffee Bean

Brisa Shiitake

* upcharge applies

BrisaShiitake

PR-515-MLA FABRICS

Model Number	PR515-SME	PR515-MLA
Size	Small/Med	Med/Large
Weight Capacity	375 LB	375 LB
Positioning	MaxiComfort	MaxiComfort
Twilight Technology	Yes	Yes
Chaise Seat	Yes	Yes
Pocketed Coil Seat	No	No
Back Type	Biscuit	Biscuit
Removable Back Option	No	No
Head Rest Cover	Yes	Yes
Arm Rest Covers	No	No
Overall Width	35.5"	39.5"
Overall Height	44"	46.5"
Width Between Arms	20"	20"
Floor to Top of Seat	20"	21.5"
Seat Depth	18.5"	20"
Seat to Top of Back	28"	31"
Seat to Top of Arm	7.5"	8.5"
Distance Required from Wall Reclined	23"	24"
Extension Reclined from Back of Seat to End of Foot Rest	43"	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes	Yes
Warranty on Electrical parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	Yes	Yes
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Customizable *	No	No
Hand Control Routing	Inside Arm	Inside Arm
Extra Pocket Option	Standard	Standard
Heat & Massage Option	No	No
Power Pillow Option	Standard	Standard
Power Adjustable Lumbar Option	Standard	Standard
Auto Flex Option	No	No
Foot Rest Extension Option	No	No
Straight Lift Option	No	No
USB Charging Port	Yes	Yes
Place of Origin	USA	USA



Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Smart Tek Our luxurious lift and

recline chairs are equipped with electronics to help assure a secure lift every time you need it.

Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.





















The World's Most Advanced **Power Positioning Lift Recliner!**

Our Patented Twilight Technology Now in a Unique 5-Motor System to Achieve Optimum Support



Motor #1

Adjusts back positioning

Controls the footrest movement and the patented movement of the rising seat box

Activates Twilight technology, a tilting movement that provides a weightless experience

NEW! Motor #4 Adjusts head support and placement.

NEW! Motor #5

Provides customizable lumbar support.

HEALTHY BENEFITS:

- Elevates Feet Above Heart Supports Natural 'S'
- Alleviates Pressure On Lungs
- Distributes Weight
- Relieves Pressure Points
- Encourages Proper Posture
- Shape of Spine
- Stylish New Track Arm Design
- Available in Luxurious **Brisa Fabrics With** Breathable Technology







